

# IMIXAgility™

Ergonomic and efficient



## IMIXAgility™



This affordable two-bucky system is made with the same attention to detail and workmanship as all of IMIX's premium systems. Agility can grow as the customer's image volume grows. Available for use with film, CR or any of the IMIX digital detectors, Agility satisfies a broad range imaging and budgetary requirements.

# IMIXAgility™ System Specifications

## TABLE

Table Design	One Column Table
Patient Load	295 kg (650 lb.)
Motorized	
Vertical Travel	54 to 85 cm (21 to 33 in.)
Longitudinal Travel	± 50 cm (± 20 in.)
Lateral Travel	±15 cm (± 6 in.)
Vertical Collision Protection	Standard
Maneuver handle	Up/Down & Table Top Release (Optional)
Foot control	Up/Down & Table Top Release (Optional)
Accessories	Supports all IMIX Intuition Table accessories

## UPRIGHT BUCKY STAND

Vertical Column Travel	147 cm. ( 58 in.) manual travel
Bucky tilt	- 20° to +90° (Optional)
Accessories	Lateral Arm Support

## FLOOR MOUNTED TUBE STAND

Column Height	237 cm (93 in.)
Focal Spot	
Height Range	40 to 184 cm (15.7 to 72 in.)
Base Length	302 cm (118 in.)
Tube Stroke	12 cm (4.7 in.)
Tube Rotation	± 135°
Pivot Base	±90° (Option)

## DETECTOR

Technology	Supports both IMIX Xwireless and Slimline technologies
------------	--

## CONTROL STATION

Image review	within 10 sec. of exposure
Display	20 in. flat panel
Image storage	Local storage for 100 images
DICOM	DICOM 3.0 complaint

## X-RAY GENERATOR

Output	65 kW high-frequency 40 kVp to 150 kVp range at 10 mA to 640 mA
Programmable	Full anatomical programming
Hand Switch	Prep/Expose hand switch

## X-RAY TUBE & COLLIMATOR

Focal Spot	0.6 to 1.2 mm
Tube Capacity	300 kHu
Target Angle	12°
Manual Collimator	Multi-blade shutter with laser line field alignment

## Option:

Rotational tubecrane movement accommodates cross table exams and flexible room configurations



\* Specifications subject to change without notice

[www.imixadr.com](http://www.imixadr.com)