

# Elite Urology Imaging Table



“Designed with  
*Your Efficiency*  
in Mind”



## **UNOBSTRUCTED C-ARM IMAGING**

- Carbon fiber tabletop
- Open design
- 4 or 5-way motorized movements

## **MULTI-PURPOSE**

- Urology
- Gynecology
- Gastroenterology
- Cardiovascular
- General Surgery
- Orthopedic
- ENT

500lbs/227kg  
Patient weight capacity



**North American Office:** 23112 Alcalde Dr. Suite A Laguna Hills, Calif. 92653  
**Worldwide Headquarters:** Anavagen 1, Vaxjo SE-325 46, Sweden

# Medstone™ Mobile Urology/GYN Imaging Table

Designed for ease of use with all mobile "C-Arm" Imaging units



Mobile urology table showing Table Top Extension and **optional** Foot and Ankle strups with disposable drain bar assembly



Patient Knee Crutches, included with system

Snap-on Patient Arm Board with pad

## 0047-BT SERIES SPECIFICATIONS

		Precision Series		
		BT-5(H)SCU	BT-4(H)SCU	BT-3ER-KOM
(1) Motorized Elevation (Lift) Travel	Minimum	26" (66.5cm)	26" (66.5cm)	26" (66.5cm)
	Maximum	38" (96cm)	38" (96cm)	38" (96cm)
(2) Motorized Trendelenburg Tilt		± 15°	± 15°	± 15°
(3) Motorized Longitudinal Table-top Travel		15¾" 40cm	15¾" 40cm	15¾" 40cm
(4) Motorized Lateral Table-top Tilt		± 15°	± 15°	NO
(5) Motorized Lateral Table-top Travel		± 4" 100mm	NO	NO
(6) Table Weights		510lbs (232kg)	490lbs (223kg)	510lbs (232kg)
(5) Shipping Weights		700lbs (318kg)	680lbs (309kg)	700lbs (318kg)
(7) Shipping Crate Size, Imperial	88" x 31" x 41"			
(8) Shipping Crate Size, Metric	230cm x 79cm x 104cm			
Carbon Fiber Table-top 24" wide x 48" long		INCLUDED	INCLUDED	INCLUDED
Carbon Fiber Table-top Extension 24" wide x 28" long		INCLUDED	INCLUDED	INCLUDED
<b>ACCESSORY KIT</b> One pair Patient Knee Crutches with Mounting Clamps Disposable Flexible Drain Bag Hoop One Starter Box of Disposable Drain Bags Hand & Foot Controls for all Movements		INCLUDED	INCLUDED	INCLUDED
Patient Weight Capacity		500lbs 227kg	500lbs 227kg	500lbs 227kg
Un-obstructed Table-top Imaging Area		58" 1473 mm	58" 1473 mm	58" 1473 mm
Radiolucent Sofcare Table top & Extension Pads		INCLUDED	INCLUDED	INCLUDED
<b>OPTIONS</b>				
Snap-on Arm Board with Mounting Clamp		OPTIONAL	OPTIONAL	OPTIONAL